

BloodCenter of Wisconsin, Inc

638 N 18 St • Milwaukee, WI 53233-2121
414.937.6250 • 800.245.3117 • Fax 414.937.6206
www.bcw.edu • labinfo@bcw.edu

Patient: **MockThromb,**

Birth Date:

Gender: Male

Alt ID:

Client: ZZ Test Organization

Hematology Genetics

More common and many rare types of inherited thrombocytopenia will be identified with this panel, including *MYH9*-related disorders, Bernard-Soulier syndrome, congenital amegakaryocytic thrombocytopenia, familial platelet disorder with predisposition to acute myelogenous leukemia, *ANKRD26*-related thrombocytopenia, *WAS*-related thrombocytopenia, gray platelet syndrome and others. Additional genes in this panel are associated with syndromes that have thrombocytopenia as a common finding among other non-hematologic features.

Assay Principle

This next-generation sequencing assay analyzes 22 genes, spanning the full coding regions plus a minimum 30bp of non-coding DNA including intron-exon junctions, in addition to approximately 200bp upstream of *ANKRD26* coding region covering all currently known intervening variants of *ANKRD26*. These targeted regions are captured by hybridization, amplified and sequenced by massively parallel sequencing. Regions will have a minimum coverage of 50x and those regions with less than 50 sequencing reads or low quality are supplemented with Sanger sequencing. All regions are covered by bi-directional analysis. Variants are identified by a customized bioinformatics pipeline, analyzed and comprehensively interpreted by our team of directors, scientists, and genetic counselors. All reported variants, including pathogenic, likely pathogenic, and variants of uncertain significance, are confirmed by Sanger sequencing.

For prenatal testing, analysis of variable number of tandem repeats (VNTR) is used to confirm results are not affected by maternal cell contamination.

Assay Sensitivity and Limitations

The analytical sensitivity of this test is >99% for single nucleotide changes and insertions and deletions of less than 20 bp. This assay does not detect large deletions or duplications (>20 bp) or deletions, duplications or variants that are outside the regions sequenced. To order the analysis of copy number variants at the exon or gene level, please refer to the aCGH Deletion/Duplication Analysis test, if available, or contact Client Services before placing your order.

This test was developed and its performance characteristics determined by BloodCenter of Wisconsin. It has not been cleared or approved by the FDA. However, the FDA has determined that such clearance is not necessary. The test has been validated in house and is used for clinical purposes. It should not be regarded as investigational or for research. Our Laboratory is certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) as qualified to perform high complexity clinical laboratory testing.

Reporting

While this assay is designed to detect germline genetic variants associated with thrombocytopenia, variants unrelated to the indication for testing, but with other clinical and/or reproductive implications, may also be detected. A comprehensive database of gene-phenotype relationships listed by gene name can be found at <http://www.omim.org>.

Hematology Genetics

Results are classified and reported in accordance with ACMG next-generation sequencing standards. Variants predicted to be pathogenic, likely pathogenic, and of uncertain significance will be reported; variants classified as likely benign or benign are typically not reported but such data are available upon request. Sequence variants are described using standard Human Genome Variation Society (HGVS) nomenclature (<http://www.hgvs.org>).

References

Balduini CL, Savoia A, Seri M. 2013. Inherited thrombocytopenias frequently diagnosed in adults. *J. Thromb. Haemost.* 11(6):1006-19.

Drachman JG. 2004. Inherited thrombocytopenia: when a low platelet count does not mean ITP. *Blood.* 103(2):390-8.

Noris P, Biino G, Pecci A, Civaschi E, Savoia A, Seri M, Melazzini F, Loffredo G, Russo G, Bozzi V, Notarangelo LD, Gresele P, Heller PG, Pujol-Moix N, Kunishima S, Cattaneo M, Bussel J, De Candia E, Cagioni C, Ramenghi U, Barozzi S, Fabris F, Balduini CL. 2014. Platelet diameters in inherited thrombocytopenias: analysis of 376 patients with all known disorders. *Blood.* 124(6):e4-e10.

Noris P, Favier R, Alessi MC, Geddis AE, Kunishima S, Heller PG, Giordano P, Niederhoffer KY, Bussel JB, Podda GM, Vianelli N, Kersseboom R, Pecci A, Gnan C, Marconi C, Auvrignon A, Cohen W, Yu JC, Iguchi A, Miller Imahiyerobo A, Boehlen F, Ghalloussi D, De Rocco D, Magini P, Civaschi E, Biino G, Seri M, Savoia A, Balduini CL. 2013. ANKRD26-related thrombocytopenia and myeloid malignancies. *Blood.* 122(11):1987-9.

Noris P, Perrotta S, Seri M, Pecci A, Gnan C, Loffredo G, Pujol-Moix N, Zecca M, Scognamiglio F, De Rocco D, Punzo F, Melazzini F, Scianguetta S, Casale M, Marconi C, Pippucci T, Amendola G, Notarangelo LD, Klersy C, Civaschi E, Balduini CL, Savoia A. 2011. Mutations in ANKRD26 are responsible for a frequent form of inherited thrombocytopenia: analysis of 78 patients from 21 families. *Blood.* 117(24):6673-80.

Pecci A, Balduini CL. 2014. Lessons in platelet production from inherited thrombocytopenias. *Br. J. Haematol.* 165(2):179-92.

Pippucci T, Savoia A, Perrotta S, Pujol-Moix N, Noris P, Castegnaro G, Pecci A, Gnan C, Punzo F, Marconi C, Gherardi S, Loffredo G, De Rocco D, Scianguetta S, Barozzi S, Magini P, Bozzi V, Dezzani L, Di Stazio M, Ferraro M, Perini G, Seri M, Balduini CL. 2011. Mutations in the 5' UTR of ANKRD26, the ankirin repeat domain 26 gene, cause an autosomal-dominant form of inherited thrombocytopenia, THC2. *Am. J. Hum. Genet.* 88(1):115-20.

Savoia A. 2016. Molecular basis of inherited thrombocytopenias: an update. *Curr. Opin. Hematol.* 23(5):486-92.

