

Custom Blood Disorder Panel

Versiti offers comprehensive genetic analysis to detect sequence variants and large deletions and duplications in any 2-10 genes across our portfolio of curated hematology genes. This analysis can be ordered as:

- Next generation sequencing (NGS) only*;
- NGS with reflex to array comparative Genomic Hybridization (aCGH) deletion/duplication analysis if sequencing does not identify clinically significant variants that fully explain the patient's phenotype;
- NGS with concurrent aCGH deletion/duplication analysis (both testing methodologies performed simultaneously); or
- Deletion/duplication analysis by aCGH only.

*Excludes *PLAU*

The custom blood disorder panel can be ordered by selecting 2-10 genes from the following curated hematology-focused genes.

Hematology Genetics Single Genes For additional information about genetic panels and more, visit: Versiti.org/HG .												
ABCG5	AP3D1	CYCS	F7	FGG	GGCX	HPS3	LMAN1	PLA2G4A	RNU4ATAC	SMAD4	TAFAZZIN	VPS13B
ABCG8	ARPC1B	DIAPH1	F8	FLI1	GIN51	HPS4	LYST	PLAU*	RUNX1	SMARCD2	TBXA2R	VPS33B
ACTB	BLOC1S3	DTNBP1	F9	FLNA	GNE	HPS5	MCFD2	PLG	SBDS	SRC	TBXAS1	VPS45
ACTN1	BLOC1S6	EFL1	F10	FYB1(FYB)	GP1BA	HPS6	MECOM	PRKACG †	SERPINA1§	SRP19	TCIRG1	VWF
ACVRL1	BTB	ELANE	F11	G6PC3	GP1BB	HRG	MPIG6B	PROC	SERPINC1	SRP54	THBD	WAS
ADAMTS13	CDC42	ENG	F13A1	GATA1	GP6	ITGA2B	MPL	PROS1	SERPIND1	SRP68	THPO	WDR1
AK2	CLPB	EPHB4	F13B	GATA2	GP9	ITGB3	MYH9	RAC2	SERPINE1	SRP72	TUBB1	WIPF1
ANKRD26	CSF3R	ETV6	FERMT3	GDF2	HAX1	JAGN1	NBEA	RASA1	SERPINF2	SRPRA	USB1	
ANO6	CXCR2	F2	FGA	GFI1	HOXA11	KDSR	NBEAL2	RASGRP2	SLC37A4	STIM1	VIPAS39	
AP3B1	CXCR4	F5	FGB	GFI1B	HPS1	KNG1	P2RY12	RBM8A	SLFN14	STXBP2	VKORC1	

aHUS/DDD Genetic Panel genes C3, C4BPA, C4BPB, CFB, CFH, CFHR1, CFHR3, CFHR4, CFHR5, CFI, DGKE, MCP are NOT available as single gene sequencing.

*PLAU available via aCGH only.

§ SERPINA1 is targeted for the Pittsburgh allele in exon 5 only.

† PRKACG NGS includes only a region to cover the p.Ile74Met variant.



Indications for testing:

Custom Blood Disorder Panel (NGS and/or aCGH), order code 4850:

The Custom Blood Disorder Panel should be considered:

- In patients with a suspected inherited disorder of coagulation factors, platelet function and/or number, thrombosis and/or neutrophils in which the clinical and/or laboratory phenotype is defined enough that the diagnostic possibilities are limited to 10 genes or less
- In patients with clinical and laboratory findings of an inherited disorder of coagulation factors, platelet function and/or number, thrombosis and/or neutrophils, when the patient's history and laboratory phenotype suggest multiple coexisting disorders
- In patients in whom a family history of an inherited hematologic disorder is reported but the specific gene or disorder is unknown, and an affected relative is unavailable for confirmation of diagnosis

Targeted Familial Variant Analysis (order code 4970):

Targeted variant analysis for clinical diagnosis, carrier identification, or prenatal diagnosis can also be performed on any gene when the pathogenic variant(s) is known in the family. If the proband was not tested at Versiti, a control sample may be needed (please call the laboratory to discuss). If the familial variant is a large deletion or duplication, analysis by aCGH for the involved gene may be required. Please call our laboratory to discuss prior to sending sample.

For clinical questions about laboratory tests and test utilization support, contact Versiti Client Services: (414) 937-6396 or 800-245-3117, Option 1, to be directed to our genetic counselors and clinical support team.

Informed Consent

It is recommended that healthcare providers obtain informed consent from the patient when genetic testing is ordered, consistent with any applicable state laws and regulations, documenting that the patient has been advised of and understands the indications for and implications of the genetic test. This panel is designed for clinical detection of germline genetic variants in genes with strong or definitive evidence for causality of hemostatic disorders. Test results may nonetheless yield genetic findings that may be unrelated to the current clinical presentation, such as risk for syndromic manifestation, predisposition or malignancy and/or may carry individual or familial implications such reproductive implications (carrier status). If needed, an informed consent form for Versiti Hematology Genetics testing can be found at <http://www.versiti.org/hg> under *forms*.

Test method:

NGS: This next-generation sequencing assay analyzes the complete coding region of 2-10 genes plus a

minimum 30bp of non-coding DNA, including intron-exon boundaries, and is compared to the build GRCh37.p13 reference sequence. These targeted regions are captured by hybridization, amplified, and sequenced by massively parallel sequencing. Regions will have a minimum coverage of 50x and those regions with less than 50 sequencing reads or low quality coverage are supplemented with Sanger sequencing. All regions are covered by bidirectional analysis. Variants are identified by a customized bioinformatics pipeline, analyzed and comprehensively interpreted by our team of practicing hematologists with expertise in non-malignant hematology and laboratory diagnostics, scientists, and genetic counselors. All reported variants, including pathogenic, likely pathogenic, and variants of uncertain significance, are confirmed by Sanger sequencing. For prenatal testing, analysis of variable number tandem repeats (VNTR) is used to confirm results are not affected by maternal cell contamination.

aCGH: The specific genes are analyzed for copy number variations due to deletion or duplication by high density gene-focused array Comparative Genomic Hybridization. Probes are approximately 60bp in length and density of coverage in exonic regions is a minimum of 4 probes per 500 bp. Genomic DNA for the samples and sex-matched references are denatured, labeled with fluorescent dye and hybridized, the array is washed and scanned, and analysis is performed for the specific genes requested.

Assay sensitivity and limitations:

NGS: The analytical sensitivity of the NGS test is >99% for single nucleotide changes and insertions and deletions of less than 20 bp. NGS analysis is not designed to detect large deletions or duplications (>20 bp), or variants that are outside the regions sequenced. Low level mosaicism will not be detected by this sequencing methodology.

aCGH: Balanced chromosomal rearrangements (i.e., translocations, inversions) or point mutations that may be the cause of the clinical phenotype cannot be detected via aCGH. Any exonic deletion or duplication smaller than 500bp may not be detected. Low level of mosaicism will not be detected by aCGH. Probe performance could be affected by multiple SNPs in a given region. Breakpoints occurring outside the targeted gene(s) will not be defined.

Clinical Sensitivity

The clinical sensitivity of the Custom Blood Disorder Panel (NGS and aCGH) is highest in patients with a clearly defined and specific clinical and laboratory phenotype for the genes selected.

Reporting of Results

Results are classified and reported in accordance with ACMG next-generation sequencing and copy number variation standards and guidelines. Sequence variants and large deletions and duplications predicted to be

pathogenic, likely pathogenic, and of uncertain significance will be reported; variants classified as likely benign or benign are typically not reported but such data are available upon request. Sequence variants are described using standard Human Genome Variation Society (HGVS) nomenclature (<http://hgvs.org>); copy number variants are described in accordance with the International System for Human Cytogenomic Nomenclature (ISCN).

Specimen Requirements

Parental/Patient/Pediatric: 3-5 mL Whole blood (EDTA tube, lavender top), 2-5 mL Bone marrow (EDTA tube, lavender top), 3-4 Buccal swabs, or $\geq 1\mu\text{g}$ of DNA at $\geq 50\text{ng}/\mu\text{L}$ of High Quality DNA.

Fetal: 7-15 mL amniotic fluid, 5-10 mg chorionic villi; back up culture of amniocytes or chorionic villi is highly recommended. Cultured: Two T25 flasks cultured amniocytes or chorionic villi (2×10^6 minimum). Maternal blood sample of 3-5 mL Whole blood (EDTA tube, lavender top) is requested for all prenatal samples for maternal cell contamination studies. For questions please contact the laboratory to discuss sample requirements.

Shipping Requirements



SHIP

Ship on an ice pack or at room temperature. Protect from freezing. Place the specimen and the requisition into plastic bags and seal. Insert into a Styrofoam container, seal and place into a sturdy cardboard box, and tape securely. Ship the package in compliance with your overnight carrier guidelines. Label with the following address:

Client Services/Diagnostic Laboratory
Versiti
638 N. 18th St
Milwaukee, WI, 53233

Required Forms



ORDER

Please complete all pages of the requisition form. Clinical history (including patient-reported ancestry, clinical diagnosis, family history, and relevant laboratory findings) is necessary for optimal interpretation of genetic test results and recommendations. Clinical and laboratory history can either be recorded on the requisition form or clinical and

laboratory reports can be submitted with the sample.

CPT Codes/Billing/Turnaround Time

Test code: 4850

CPT codes: Visit [Versiti.org/test-catalog](https://www.versiti.org/test-catalog) for full list of recommended CPT codes for each gene in the panel.

Turnaround time: NGS only, aCGH only, or NGS and aCGH concurrently: 21 days

NGS reflex to aCGH: 21 days (if NGS only, aCGH not needed) or 42 days (with reflex to aCGH)

The CPT codes provided are subject to change as more information becomes available. CPT codes are provided only as guidance to assist clients with billing.

For additional information related to shipping, billing or pricing, please contact Versiti Client Services: (414) 937-6396 or 800-245-3117, Option 1, or LabInfo@versiti.org

References

Additional references for inherited coagulation, platelet, neutrophil disorders and thrombosis can be found under their specific panel test descriptions

Variant interpretation references

1. Bean LJH, Funke B, Carlston CM, et al. Diagnostic gene sequencing panels: from design to report-a technical standard of the American College of Medical Genetics and Genomics (ACMG). *Genet Med.* 2020;22(3):453-461. doi:10.1038/s41436-019-0666-z
2. Rehm HL, Bale SJ, Bayrak-Toydemir P, et al. ACMG clinical laboratory standards for next-generation sequencing. *Genet Med.* 2013;15(9):733-747. doi:10.1038/gim.2013.92
3. Richards S, Aziz N, Bale S, et al. Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. *Genet Med.* 2015;17(5):405-424. doi:10.1038/gim.2015.30
4. Riggs ER, Andersen EF, Cherry AM, et al. Technical standards for the interpretation and reporting of constitutional copy-number variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics (ACMG) and the Clinical Genome Resource (ClinGen). *Genet Med.* 2020;22(2):245-257. doi:10.1038/s41436-019-0686-8.

