

ALERT

Known Historical Mismatch Sample Arriving

(Pre-Testing)

Facility / Client: _____

Notified by: _____

Sample ID #: _____

Collection Date: _____

Historical Result: A Pos B Pos O Pos AB Pos A Neg B Neg O Neg AB Neg

Historical Phenotype: C+ Cneg c+ cneg E+ Eneg e+ eneg

K+ K neg

Sample Arriving: _____

Email completed form to TLLR@indianablood.org

For Versiti DTL Use Only:

Lot Release: place the sample on “Hold”.

After receiving current test results, verify the historical result matches the current result.

- If results are discrepant, forward the tube and testing results to IRL.
- If results match, transfer the current test results to the client.